

Mark X™

User Manual





Table of Contents

Mark X™ User Manual

Important Safety Information	i
Specifications	iii
Chapter 1: Mark X™ System	1
The Mark X™ Color Panel Options	2
The Mark X™ Labels	3
The Mark X™ Features	4
The Mark X™ Important Part Numbers	5
The Mark X™ Key Connections Overview	7
Connecting the Detector Cable f/Xmaru 1210 PCB	8
Connecting the Generator	9
Powering Up the Mark X™ System	10
Chapter 2: Opening Mark X™ Software	11
Chapter 3: Adding Users	13
Chapter 4: Patient Information	14
Searching and Adding Patients	14
Editing an Existing Patient's Information	16
Creating a New Study for an Existing Patient	16
Chapter 5: Creating Radiographs	17
Choosing Anatomy and Positioning	17
Acquiring Radiographs	19
Cropping a Radiograph	21
Chapter 6: Protocols	23
Using Protocols	23
Adding New Protocols	26
Editing Existing Protocols	30
Renaming Protocols	31
Deleting Protocols	33



Table of Contents

Mark X™ User Manual *(continued)*

Chapter 7: Managing Studies	35
Adding New Images to an Existing Study.....	35
Reviewing a Study	36
Chapter 7: Managing Studies	37
Chapter 8: Exporting Images to a Folder	39
Chapter 9: Transferring Studies to the Storage Station.....	41
Transferring Individual Studies to the Storage Station.....	43
Transferring DICOM Batch to the Storage Station	44
Confirm Storage of a Study	45
Chapter 10: Adding and Editing a DICOM Location.....	46
Appendix A: Digital Radiography Technique Charts	49
Portable 80/15	49
Portable 100/30.....	52
Gantry	54
Appendix B: Mark X™ System Field Replacement Guide	56
Appendix C: Calibration Procedure	60
Appendix D: Technical Support.....	62

Important Safety Information

Note



*Read all the instructions before connecting or operating the component
Keep this manual so you can refer to these safety instructions.*

Warning



*There are no user serviceable parts inside.
Refer all servicing to qualified service personnel.*

To reduce the risk of fire or electric shock, do not expose the unit to moisture or water. Do not allow foreign objects to get into the enclosure. If the unit is exposed to moisture, or a foreign object gets into the enclosure, immediately disconnect the power cord from the wall. Take the unit to a qualified service person for inspection and necessary repairs.

Caution



Heed all warnings and safety information in these instructions and on the product itself. Follow all operating instructions.

*Clean the enclosure only with a dry cloth or a vacuum cleaner.
Clean the DR plate with a damp cloth. Use no chemicals or alcohol.*



IPX0 per IEC 60529

Type B per IEC 60601-1

You must allow 10 cm or 4 inches of unobstructed clearance around the unit. Do not place the unit on a bed, sofa, rug, or similar surface that could block the ventilation slots. If the component is placed in a bookcase or cabinet, there must be ventilation of the cabinet to allow proper cooling.

Keep the component away from radiators, heat registers, stoves, or any other appliance that produces heat.

Keep the component away from flammable materials.

The unit must be connected to a power supply only of the type and voltage specified on the side panel of the unit

Important Safety Information *(continued)*



Connect the component to the power outlet only with the supplied power supply cable or an exact equivalent. Do not modify the supplied cable in any way. Do not attempt to defeat grounding and/or polarization provisions. Do not use extension cords.

Do not route the power cord where it will be crushed, pinched, bent at severe angles, exposed to heat, or damaged in any way. Pay particular attention to the power cord at the plug and where it exits the back of the unit.

The power cord should be unplugged from the wall outlet if the unit is to be left unused for a long period of time. Immediately stop using the component and have it inspected and/or serviced by a qualified service agency if:

- The power supply cord or plug has been damaged.*
- Objects have fallen or liquid has been spilled into the unit.*
- The unit has been exposed to rain.*
- The unit shows signs of improper operation*
- The unit has been dropped or damaged in any way*

Place the unit on a fixed, level surface strong enough to support its weight. Do not place it on a moveable cart that could tip over.

System will become unresponsive when exposed to Electrostatic Discharge (ESD) or Electrical Fast Transient (EFT) conditions. This is the normal response to these environmental conditions. Please restart the system to return to normal operation

USB ports are for keyboard and mouse only. User is responsible for evaluating the system with any devices connected to this port for compliance to IEC 60601-1-1.



Ethernet port is not to be used during acquisition. User is responsible for evaluating the system with any devices connected to this port for compliance to IEC60601-1-1.

Specifications

*Operation Environment: +15°C to 35°C,
30% to 75% RH, 70 to 106kPa (pressure)*

*Storage / Transport : -10° C to 55° C, 15% to 80% RH,
70 to 106kPa (pressure)*

*Power Requirements:
100 - 240 VAC 5.6A 50/60Hz
Weight (net) 47 lbs*




CAUTION

Risk of electric shock.
Do not open.


Caution: To reduce the risk of electric shock, do not remove cover. No user-serviceable parts inside. Refer servicing to qualified service personnel.

Applicable for USA, Canada or where approved for the usage.

Caution: To prevent electric shock, match wide blade of plug to wide slot. Insert fully.



This symbol is to alert the user to important operating and maintenance (service) instructions in this manual and literature accompanying the product.



This symbol is to alert the user to the presence of uninsulated dangerous voltages inside the product's enclosure that may constitute a risk of electric shock.

Chapter 1: Mark X™ System

The Mark X™ system links seamlessly with the X-ray generator to provide the best image quality in a small, lightweight, self-contained system. The Mark X™ portable Acquisition Stations are intended to be used by Veterinarians for intermittent mobile diagnostic applications. They will connect the Mark X Device to a Xmaru 1210PCB Amorphous Silicon Digital X-Ray Detector.



Chapter 1: Mark X™ System *(continued)*

The Mark X™ Color Panel Options



30-901 Red



30-900 Purple



30-904 Silver



30-902 Blue



30-899 Gun Metal



30-898 Green

Chapter 1: Mark X™ System *(continued)*

The Mark X™ Labels



Chapter 1: Mark X™ System *(continued)*

The Mark X™ Features



Chapter 1: Mark X™ System *(continued)*

The Mark X™ Important Part Numbers



94-828 Mark X™



34-828 Mark X™,
Electronics Box Assembly



34-827 Mark X™,
A-Frame Assembly

Chapter 1: Mark X™ System

The Mark X™ Important Part Numbers (*continued*)



10-121 Min-X-Ray Sync Cable for Mark X™
Syncs the generator's prep and expose cycle with the computer.



90-403 10' Ethernet Patch Cable
Connects the computer to the network to allow for easy image transfers from the Mark X™ to PACS.



91-408 Xmaru 1210PCB Detector
Captures digital radiographs.



10-124 22' Detector Cable f/Xmaru 1210PCB
Provides power to the plate and syncs the plate and Mark X™ System.



10-171 Equine 20' IEC-NEMA Power Extension Cord
Connects the Mark X™ to any standard U.S. 110V power outlet (Nema 5-15).

Chapter 1: Mark X™ System *(continued)*


The Mark X™ Key Connections Overview



Chapter 1: Mark X™ System *(continued)*

Connecting the Detector Cable f/Xmaru 1210 PCB

1. Remove the Plate and the Plate Cable.
2. Connect the 22' Detector cable f/Xmaru 1210PCB (10-124) by aligning the red dot and notch on the cable with the red dot and notch on the receptacle.
3. Connect plate cable to the plate if it is not already connected.



Chapter 1: Mark X™ System *(continued)*

Connecting the Generator

The sync cable should stay connected to the generator and should be stored in the generator box.

1. Plug the Min-X-ray Sync Cable (10-121) into the Mark X™. Align by red dots on each end of the connector.
2. Plug the hand switch end of the Min-X-ray Sync Cable (10-121) into the generator hand switch port.



Chapter 1: Mark X™ System *(continued)*


Powering Up the Mark X™ System

1. Plug the female end of the power cable into the power cable slot on the Mark X box. Plug the other end into a grounded power outlet (see page ii Power Requirements).
2. With the cables connected, turn on the Mark X™ by pressing the power button on the top of the unit.



Chapter 2: Opening Mark X™ Software

Open the Mark X™ software to begin the process of taking a radiograph.

1. Double-tap the TruDR™ icon  on the desktop. This will load the TruDR™ software for acquiring digital radiographs.
2. Enter the user name and password when the login screen appears.

The default login is:

User Name: admin

Password: admin



The software will load directly into the Patient Screen if the “Remember Me” box is checked.






Chapter 3: Adding Users

Add veterinarians and technicians to TruDR™.


1. Touch the **Manage** button  in the Acquisition Screen.



To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 1.800.268.5354.



2. Select the **Users** menu from the left side of the screen to load the Users Screen.



Chapter 3: Adding Users *(continued)*

3. Touch the **Add User** button in the Acquisition Screen.
4. Enter the User's name in the **First Name** and **Last Name** fields.
5. Create a unique user name and password in the **User Name** and Password fields.



Suggested user name: First initial of your first name plus last name (e.g. jsmith for John Smith). Passwords can be left blank for ease of use.

6. Select a Group or Groups for each user by Touching **Tech** or **Vet**.

Each Group allows different permissions. You may select multiple check boxes per user:


Tech: User is added to the Tech pull down list has no access to the Manage features.

Vet: User is added to the Vet pull down list, has no access to the Manage features.

Admin: Has full access to the Manage features.

Tech Support: Should only be used by Sound-Eklin™ Tech Support

7. Touch the Save button  to add the User.



User ID	First Name	Last Name	Email	Group
3	Admin	First Name	Admin@soundeklin.com	Admin
7	Jackie	Chiles	jchiles@soundeklin.com	Vet, Admin
5	Joe	David	jedavid@soundeklin.com	Tech
6	Henry	David	hmdavid@soundeklin.com	Vet, Admin
4	TechSupport	First Name	TechSupport@soundeklin.com	TechSupport

Selected User Info

User ID: New User

First Name: User Name:

Last Name: Password:

Email:

Groups: ☐ Tech ☐ Vet ☐ Admin ☐ TechSupport

8. Repeat steps 3 through 7 to add additional Users.



Please note: There must be at least one TECH active in the user menu.

Chapter 4: Patient Information


Searching and Adding Patients

Search for existing patients or add new patients easily by using the filtered search fields.

1. Search for a patient in the Patient Screen by entering the **Patient ID**, the **Patient Name**, or the **Owner Last Name**. The patient list will automatically filter as you type.



Patient list is displayed based on the User Selected Default.

2. If you do not see the patient you are searching for, Touch the Add  button or press the **Enter** key on the keyboard to bring up the Add Patient Screen.



3. Enter the patient's information. **Patient Name** and **Last Name** are required fields.



When entering patient name and last name, DO NOT use special characters such as: % \$ &. Special characters and symbols can cause problems in the patient name fields.

Chapter 4: Patient Information

Searching and Adding Patients (*continued*)

4. Touch one of the Save buttons.



Save

Saves the current patient and returns you to the Patient Screen.



Save+ Add

Saves the current patient and returns to the Add Patient Screen.



Save+ Acquire

Saves the current patient then opens the Anatomy/Protocol Screen.



Save+ Protocol

Saves the patient and brings you to the Protocol Screen to select the desired series.

The screenshot shows a tablet interface for patient information. At the top, there are several icons: 'patient', 'save', 'save+add', 'save+acquire', 'save+protocol', and 'keyboard'. Below these icons is a list of patients, with 'Jones' selected. The main form is divided into two sections: 'Patient Information' and 'Owner Information'. The 'Patient Information' section includes fields for ID, F Name, L Name, Sec. ID, Species, Breed, Sex, Weight, and DOB. The 'Owner Information' section includes fields for Last Name, First Name, Institution, Address, Address2, Address3, and City. There are also buttons for 'Save', 'Save+ Add', 'Save+ Acquire', 'Save+ Protocol', and 'Clear Owner'. A keyboard is visible at the bottom of the screen.

Chapter 4: Patient Information *(continued)*

Editing an Existing Patient's Information

Edit an existing patient's information for future DICOM tagging.



Previous exam information cannot be changed

1. Select a patient to edit from the Patient Screen.
2. Touch the **Edit** button.



3. Edit the information as needed.



If the patient is from an owner previously entered, owner can be selected from owner list on the left hand side.

4. Touch the appropriate **Save** button (*see Page 17*).

Creating a New Study for an Existing Patient

Create a new study for a patient previously entered in TruDR™.

1. Search for a patient in the Patient Screen by entering the Patient ID, the Patient Name, or the Owner Last Name. The Patient List will automatically filter as you type.
2. Select a patient by Touching on the name in the Patient List.
3. Touch the **New** button to create a new study.




4. Select the anatomy and the view from the Anatomy Screen (*see Page 19*).

Chapter 5: Creating Radiographs

Choosing Anatomy and Positioning


Choose the anatomy and positioning to automatically enhance and position the radiograph.

1. Select the **Tech** and ordering **Vet** from the drop-down menus.



A screenshot of a software interface showing two drop-down menus. The first menu is labeled 'Tech:' and has 'Davola, Joe' selected. The second menu is labeled 'Vet:' and has 'Bania, Kenny' selected. Both menus have a small blue arrow icon on the right side.


2. Touch the Mark X™ touchscreen and select the appropriate anatomy.



Chapter 5: Creating Radiographs

Choosing Anatomy and Positioning (*continued*)

3. Touch on the appropriate view button for the exposure.
This will bring you to the acquisition screen.



Chapter 5: Creating Radiographs *(continued)*

Acquiring Radiographs

1. Set the X-ray generator according to the Sound-Eklin™ Technique Chart using the appropriate kvp and mAs (depending on the X-ray generator).



Due to variations between X-ray generators, optimal techniques may vary somewhat from the values printed on the Sound-Eklin™ Technique Chart.

Users should be fully acquainted with State and local regulations governing radiation protection and the use of diagnostic X-ray equipment as well as the manual for the X-Ray Generator being used.

2. Align the generator 22 inches from the plate using the light beam to fill the entire capture area defined on the plate.



If the generator is aligned perpendicular, the light will form a rectangle with even sides filling the capture area.


The image capture area is highlighted by a white box with the dotted cross-hairs in the middle.

3. From the Acquisition Screen press the hand switch halfway down to prep and all the way down to expose the X-ray machine.



Proper Two-Stage Hand Switch Technique


1. To Prep, push the hand switch halfway down.
2. Wait until the Generator is ready to fire (i.e. prep light on the X-ray machine is illuminated, etc.).
3. To expose, push the hand switch completely down .



Chapter 5: Creating Radiographs

Acquiring Radiographs (*continued*)

- The image will appear in less than 3 seconds. Use the Use Left, Right, Flip, Reverse, ROI or Crop buttons to make orientation adjustments before shooting the next shot.





'Window' and level with finger in up and down, left or right motion.

Crop the images as needed (see Pages 23-24).

If you need to reshoot your X-ray, press the hand switch and fire again.





- Touch the **Anatomy** button  to change views.
- Select the next view from the Views Screen. Touch the **Anatomy** button from the Views Screen if another anatomy is needed.
- Touch the **Review** button  to review. Touch the Patient button to return to the Patient Screen.

Chapter 5: Creating Radiographs *(continued)*

Cropping a Radiograph

Remove unwanted portions of the radiograph.

1. Make adjustments before shooting the next radiograph. To do so, Touch the Action button.  Touch the Crop button  in the Acquisition Screen.



*The new image will appear lighter after Touching the crop button.
Based on the crop selection the image will be windowed and leveled again.*



Chapter 5: Creating Radiographs

Cropping a Radiograph (*continued*)

2. Frame the desired image dragging from upper left corner to lower right hand corner of the touchscreen with finger.



3. The image will enhance once it is cropped. The image will perform an ROI based on crop size.
4. Re-crop the image if necessary before taking the next radiograph by repeating the steps.

Chapter 6: Protocols

Using Protocols

Shoot a series of X-rays in a preset order which automatically advances to the next shot after each shot taken.

1. From the Patient Screen, select the **Patient** and Touch **Protocol**.



Or, from the Patient Screen, Touch **New** and Touch **Protocol**.




OR, from the Add Patient or Edit Patient screen, Touch **Save + Protocol**.
(to Add a Patient see page 16).




Chapter 6: Protocols

Using Protocols (*continued*)

2. Select a Protocol button to begin the protocol series.





3. The Acquisition Screen now has a Protocol Information bar showing the Protocol shot number. The current Shot Information is located in the upper left portion of the screen.



Chapter 6: Protocols

Using Protocols (*continued*)

4. Shoot the first image. The software will automatically advance to the next shot in the protocol series. The shot number in the Protocol Information Bar keeps track of the shot sequence, while the current Shot Information appears in the upper left.



5. Use the Protocol Information bar to repeat a shot, skip a shot, restart the protocol, move to the last shot in the protocol, or stop the protocol.



Move to the first shot in the protocol series



Repeat the previous shot




Skip to the next shot



Skip to the last shot in the protocol series





Stop the protocol

6. Touch the **Review** button to open the Review screen (*see Page 25*).
7. Touch the **Patient** button  to return to the Patient Screen.

Chapter 6: Protocols (continued)


Adding New Protocols

Create new protocols for your individual studies.

1. Touch the Actions button  to access the Manage button  in the Patient Screen.





To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 1.800.268.5354.



Chapter 6: Protocols

Adding New Protocols (*continued*)

2. Touch on the **Plus Sign**  next to the Protocols List to show the existing protocols.
3. Select the **Protocols** menu from the left side of the screen to load the Protocol Screen. Navigate to Species > Equine > 1210S CIE > Protocol.




4. Enter the name of your new protocol in the **Protocol Name** field.
5. Touch the **Save** button. 




Chapter 6: Protocols

Adding New Protocols (*continued*)

6. Filter the shots as necessary by selecting the **Filter by Laterality** drop down menu.




7. Select the first shot from the Shots for Species list.








8. Touch the **Link** button <<Link to add the shot to the Shots Linked to Protocol box.



Chapter 6: Protocols

Adding New Protocols (*continued*)




9. Repeat step 6 through 8 as necessary.
10. Touch the **Unlink** button  to remove any shots from the **Shots Linked to Protocol** box.
11. Select a shot in the **Shots Linked to Protocol** box and change the order using the **Move Up**, **Move Down**, **Move First**, and **Move Last** buttons.

	Moves your protocol shot up one shot in the series.
	Moves your protocol shot down one shot in the series.
	Moves your protocol shot from its location to the first shot in the series.
	Moves your protocol shot from its location to the last shot in the series.

12. Touch the **Save** button  to save your protocols.
13. Touch the **Patient** button  to return to the Patient Screen.





If you have forgotten to save your protocols, a Save dialogue box will appear. Touch the Yes button to save your changes and return to the Patient Screen.



Chapter 6: Protocols (continued)


Editing Existing Protocols


Edit current Protocols to customize them to your shot series.


1. Touch the Actions button  to access the Manage button  in the Patient Screen.



To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 1.800.268.5354.




2. Touch on the Plus Sign  next to the Protocols list to show the existing protocols.
3. Select a Protocol to edit from the left side menu, make changes and Touch Save.



Chapter 6: Protocols

Editing Existing Protocols (*continued*)



- Continue with Step 8 on Page 34.



Closing the software and re-opening is required to save protocol database changes.


Renaming Protocols

Rename current protocols.

- Touch the Actions button  to access the Manage button  in the Patient Screen.





To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 1.800.268.5354.




Chapter 6: Protocols

Renaming Protocols (*continued*)

2. Touch on the **Plus Sign**  next to the Protocols List to show the existing protocols.
3. Select the **Protocol** to rename from the left side menu. Touch the **Rename** button.



4. Rename the protocol in the text box and Touch the **OK** button when finished.





5. Restart the software by touching Exit from the Main Patient Screen



Chapter 6: Protocols (continued)


Deleting Protocols


Delete current protocols.


1. Touch the Actions button  to access the Manage button  in the Patient Screen.



To access the Manage screen, the user needs to have Administrative rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 1.800.268.5354.



2. Touch on the Plus Sign  next to the Protocols List to show the existing protocols.
3. Select the Protocol to delete from the left side menu Touch the Delete button.



Chapter 6: Protocols

Deleting Protocols (*continued*)

4. Touch the **OK** button  to permanently remove the protocol.





5. Restart the software by Touching the close button in the upper right corner.

Chapter 7: Managing Studies


Adding New Images to an Existing Study

Add new images to an existing patient's study.

1. Search for a patient in the Patient Screen by entering the **Patient ID**, the **Patient Name**, or the **Owner Last Name**. The Patient List will automatically filter as you type.
2. Touch the Acquire button. 




This will take you to the Views Screen using the anatomy you selected from the original study. This is shown in the description field.




Chapter 7: Managing Studies (continued)

Reviewing a Study

Change the layout or delete any unwanted images before the study is archived on the Storage Station.

1. Select a patient by Touching on the name in the **Patient List**.
2. Touch on the **Study** from the list on the right.
3. Touch the **Review** button. 



4. Select a thumbnail in the left-hand pane to load the image into the review pane on the right. The selected image is highlighted by a yellow box.



Chapter 7: Managing Studies

Reviewing a Study (*continued*)

5. Review the images by using the following techniques:

Zoom: Once the Zoom button is selected tap the screen on area you would like to Zoom in on. To Zoom out tap revert button.

Pan: Tap the Pan button, Touch and drag the image with your finger.

Window Level: Tap Window Level button, Move up and down with finger on screen to perform Level. Dragging finger Left to Right is Window.

6. Use the toolbar for advanced review features.



Patient: Returns to the Patient Screen.



New: Opens the Anatomy/Protocol Screen for creating a new study or protocol.



Acquire: Opens the Views Screen using the anatomy you selected from the original study.



Anatomy: Opens the Anatomy/Protocol Screen to select a new view to add into the same study.



Save: Saves the image with any changes that have altered the original image.



Warning message will appear to choose “Yes” to overwrite (replace) the image or “No” to create a new (additional) image.



Delete: Removes any unwanted images.



Deleting an image permanently removes it from the computer.

Chapter 7: Managing Studies

Reviewing a Study (*continued*)



Revert: Reverts to the original image.



Left: Rotates the image 90° counterclockwise.



Right: Rotates the image 90° clockwise.



Flip: Flips the image vertically.



Reverse: Flips the image horizontally.



Magnify: Activates the magnifier which will zoom on the selected area.

Actions: Shows a drop-down menu with additional review features.



Win Level: Selects the window level feature (also accessible by using the right mouse button).



Invert: Inverts blacks and whites, replicating a fluoroscopy.



Pan: Pans the image (also accessible by using the left mouse button).



Crop: Trims an image without applying enhancements.

In the Review Screen, cropping does not re-enhance the image based on the captured area. It only trims the image without re-enhancing the image.



Overlay: Opens menu for DICOM text display.



Hide: Hides all DICOM information.



Summary: Displays patient and shot information.



Detail: Displays patient, shot, and additional detailed information.





ROI: Draw a box from left to right to window and level based on the Region of Interest.

7. Touch the **Patient** button  to return to the Patient Screen.

Chapter 8: Exporting Images to a Folder

Export images from a study to a folder to save on a computer, thumb drive, etc.

1. Search for a patient in the Patient Screen by entering the **Patient ID**, the **Patient Name**, or the **Owner Last Name**. The Patient List will automatically filter as you type.
2. Select a patient and a study. The patient and study are highlighted in a blue box.
3. Touch the **Folder** button. 



4. Select either the DICOM or JPEG image format from the Export As: drop-down menu.



5. When selecting the DICOM option from Step 4, Touch the Add DICOM viewer check box only if the computer does NOT have a full DICOM viewer already installed.



Annotation Burning
Anonymize Images
Repository Format
Create DICOM DIR

Imprints annotations made to saved images
Strips all client and clinic information from the images
Places images in a repository folder format
Adds a DICOM directory to quickly import all images into a PACS system.

Chapter 8: Exporting Images to a Folder *(continued)*

- When selecting the JPEG format from Step 4, select the Overlay Burning options from the drop-down menu:



The default option is set to None:

- None:** Does not show overlays on the JPEG image
- Detail:** Shows a detailed DICOM overlay on the JPEG image (shows Image Number, Manufacturers Model Number, Study ID, Image Anatomy, View Position, Image Laterality Extended, Columns X Rows, Technique, Institution Name, Physician of Record, Patients Name, Patient ID, Patient Breed Description, Patients Birth Date, Patients Sex, Acquisition Number, Acquisition Date, Acquisition Time.)
- Summary:** Shows a summary DICOM overlay on the JPEG image (shows Image Number, Image Anatomy, View Position, Image Laterality Extended, Institution Name, Physician of Record, Patients Name, Patient Breed Description, Patients Birth Date, Patients Sex, Acquisition Date, Acquisition Time).

- Touch the Select Button



- Select a location to save the images.



Images are saved to a selected location in a folder with the animals name and the clients last name (e.g. Trigger-Smith).



Chapter 9:

Transferring Studies to the Storage Station

Transfer the studies from the Acquisition Station to the Storage/Review Station.



You must have your Mark X™ and your Storage/Review Station networked. Please contact your local IT Person for any networking needs and SoundEklın tech support for configuring DICOM Server for the Mark X.



Unit cannot be connected to a network while in use. In addition you must remove network plug



1. Connect a Cat-5 Internet cable into the Internet port on the computer.




2. To bring up the Store Files to DICOM Server screen, Touch the Server Button  on the Patient Screen.



Chapter 9: Transferring Studies to the Storage Station *(continued)*

3. Touch the Start Store button to begin the DICOM image transfer to the Storage Station.



This will store all the studies that were taken since your last transfer and clear the queue. If any studies are missing, you can transfer them individually (see the following section).


4. Status box displays transfer results. When the status box states “completed”, your images have successfully transferred to the DICOM server. Touch the Red close button in the upper right corner to close the screen.




Chapter 9: Transferring Studies to the Storage Station *(continued)*

Transferring Individual Studies to the Storage Station

Transfer the studies from the Acquisition Station to the Storage/Review Station.

1. Search for a patient's study by entering the Patient ID, the Patient Name, or the Owner Last Name. The Patient List will automatically filter as you type.
2. Select a patient by Touching on the patient name in the Patient List field.
3. Touch the Server button. 



4. Touch the Start Store button  to begin DICOM image transfer.



Chapter 9: Transferring Studies to the Storage Station


Transferring Individual Studies to the Storage Station (*continued*)

5. Status box displays transfer results. When the status box states “completed”, your images have successfully transferred to the DICOM server. Touch the Red close button in the upper right corner to close the screen.



Transferring DICOM Batch to the Storage Station

1. If you wish to send a group of selected studies from a specific timeframe, enter the start and end date of these studies in the display at the top of the Batch Send Options window.. The Dicom image transfer will send theses studies listed in the queue all at once.



2. Touch the Start Store button **Start Store** to begin DICOM image transfer.

Chapter 9: Transferring Studies to the Storage Station *(continued)*

Confirm Storage of a Study

Confirm the study successfully sent to the Storage/Review Station.


1. On the Patient Screen, the Archive field will show a server name (e.g. SOUND-SERV1).
2. Hold the mouse over the server name to show a pop-up window with the transfer data.




The pop-up box displays how many images have been transferred and the date and time of the transfer. The information is color coded for easy reference. If only one study appears, it will be highlighted in blue and the color coding will not appear:

Green Box— Image transfer complete

Red Box — Image transfer failed



Notice the Secondary Transfer field. Secondary Transfer documents the number of studies sent to another PACS system besides the main archive, ie Sound-Serv2.



Chapter 10:


Adding and Editing a DICOM Location

Add DICOM server locations for sending images to a Storage Station, Radiologist, etc.

1. Touch the **Manage** button  in the Acquisition Screen.



To access the Manage Screen, the user needs to have Admin rights. Changing items in the Manage Screen other than “Users” and “Protocols” may cause your system to function improperly. If you are concerned about making changes in the Manage Screen, please contact Technical Support at 1.800.268.5354.




2. Select the **Configuration** menu from the left side of the screen to load the Users screen.



Chapter 10: Adding and Editing a DICOM Location *(continued)*

3. Select the **DICOM** button from the left to bring up the DICOM Server Screen.



4. Touch the Add button **Add** to add a new DICOM server. Or if editing select the existing DICOM server from the DICOM Servers list on the Right.
5. Enter the **AE Title**, **IP Address**, and **Port** number for your DICOM server. OR if editing enter the corrected **AE Title**, **IP Address**, and **Port** number for your DICOM server.



Contact Sound-Eklin™ Technical Support or the location you are sending the DICOM images for the appropriate DICOM provider information.

6. Select additional information by choosing Yes from the drop-down menus.



Only select the following options if you are directed to do so by a Technical Support representative:

Auto Send Server: Automatically sends images to the selected DICOM location.

Allow Worklist: Designates the DICOM location as a worklist modality.


Telemedicine Server: Designates the DICOM location as a telemedicine server.


Archive Server: Designates the DICOM location as an archive server.

7. Touch the **Save** button **Save** to save the current information.

Chapter 10: Adding and Editing a DICOM Location *(continued)*

8. Touch the **Verify** button  to verify that you can connect to the DICOM Server.



9. A smiley face  will appear when you have verified the connection to the DICOM Server.



*If a frowning face appears, please contact Sound-Eklin™
Technical Support (1.800.268.5354).*

10. Touch the **Patient** button  to return to the Patient Screen.

Appendix A:

Digital Radiography Technique Charts

Portable 80/15

Study	View	SID in/cm	kVP	Time in Seconds
Navicular				
	65 Degree Navicular	24/61	65	0.2
	65 Degree P-3	24/61	80	0.05
	DP	24/61	80	0.1
	Lateral	24/61	80	0.1
	Flexor Skyline	18/46	80	0.1
If Indicated	Nav. Obliques	24/61	80	0.1
If Indicated	65 P-3 Obl.	24/61	80	0.05
P-3	Lateral	24/61	80	0.1
	65 Degree P-3	24/61	80	0.05
If Indicated	65 Degree Obliques	24/61	80	0.05
Pastern				
	DP	24/61	80	0.1
	Lateral	24/61	80	0.1
	Medial Oblique	24/61	80	0.1
	Lateral Oblique	25/64	80	0.1
Fetlock				
	DP	24/61	80	0.1
	Medial Oblique	24/61	80	0.1
	Lateral	24/61	80	0.1
	Lateral Oblique	24/61	80	0.1
	Flexed Lateral	24/61	80	0.1
If Indicated	Special Palmar	24/61	80	0.1
Metacarpus/Metatarsus				
	DP	24/61	80	0.07
	Lateral	24/61	80	0.07
	Medial Oblique	24/61	80	0.07
	Lateral Oblique	24/61	80	0.07

Appendix A: Digital Radiography Technique Charts

Portable 80/15 (continued)

Study	View	SID in/cm	kVp	Time in Seconds
Carpus				
	DP	24/61	80	0.1
	Medial Oblique	24/61	80	0.1
	Lateral	24/61	80	0.1
	Lateral Oblique	24/61	80	0.1
	Flexed Lateral	24/61	80	0.1
If Indicated	Distal Rad. Sky	18/46	80	0.1
	Distal Row Sky	18/46	80	0.1
If Indicated	Prox. Row Sky	18/46	80	0.1
Tarsus				
	DP	24/61	80	0.15
	Medial Oblique	24/61	80	0.1
Angle Down 5	Lateral	24/61	80	0.1
Angle Up 5	Lateral Oblique	24/61	80	0.1
If Indicated	Calcaneal Skyline	18/46	80	0.1
Skull (TMJ & Bullae)				
	DV	24/61	80	0.25
	Lateral	24/61	80	0.25
	Right Oblique	24/61	80	0.25
	Left Oblique	24/61	80	0.25
Sinus/Dental				
	DV	24/61	80	0.2
	Lateral	24/61	80	0.2
	Right Oblique	24/61	80	0.2
	Left Oblique	24/61	80	0.2
Cervical Spine				
	Cranial (Lateral)	24/61	80	0.2
	Middle (Lateral)	24/61	80	0.25
	Caudal (Lateral)	24/61	80	0.3

Appendix A: Digital Radiography Technique Charts

Portable 80/15 (*continued*)

Study	View	SID in/cm	kVp	Time in Seconds
Pharynx	Lateral	24/61	80	0.25
Shoulder	Lateral	24/61	80	0.4
Elbow	Cran->Caud	24/61	80	0.3
	Medio-Lateral	24/61	80	0.2
Radius				
	Cran->Caud	24/61	80	0.2
	Lateral	24/61	80	0.2
Stifle				
	Caud->Cran	24/61	80	0.3
	Lateral	24/61	80	0.2
	Lateral Oblique	24/61	80	0.2
Tibia				
	Caud->Cran	24/61	80	0.2
	Lateral	24/61	80	0.2
Withers	Lateral	24/61	80	0.1

Appendix A: Digital Radiography Technique Charts *(continued)*

Portable 100/30

Study	View	SID in/cm	kVp	mAs
Navicular				
	65 Degree Navicular	30/76	80	3.2
	65 Degree P-3	30/76	80	0.5
	DP	30/76	80	1.7
	Lateral	30/76	80	1.7
	Flexor Skyline	18/46	80	1.7
If Indicated	Nav. Obliques	30/76	80	1.7
If Indicated	65 P-3 Obl.	30/76	80	0.5
P-3				
	Lateral	30/76	80	1.7
	65 Degree P-3	30/76	80	0.5
If Indicated	65 Degree Obliques	30/76	80	0.5
Pastern				
	DP	30/76	80	1.7
	Lateral	30/76	80	1.7
	Medial Oblique	30/76	80	1.7
	Lateral Oblique	30/76	80	1.7
Fetlock				
	DP	30/76	80	1.7
	Medial Oblique	30/76	80	1.7
	Lateral	30/76	80	1.7
	Lateral Oblique	30/76	80	1.7
	Flexed Lateral	30/76	80	1.7
If Indicated	Special Palmar	30/76	90	1.7
Metacarpus/Metatarsus				
	DP	30/76	80	1.7
	Lateral	30/76	80	1.7
	Medial Oblique	30/76	80	1.7
	Lateral Oblique	30/76	80	1.7
Carpus				
	DP	30/76	80	1.7
	Medial Oblique	30/76	80	1.7
	Lateral	30/76	80	1.7
	Lateral Oblique	30/76	80	1.7
	Flexed Lateral	30/76	80	1.7
If Indicated	Distal Rad. Sky	18/46	80	1.7
	Distal Row Sky	18/46	80	1.7
If Indicated	Prox. Row Sky	18/46	80	1.7

Appendix A: Digital Radiography Technique Charts

Portable 100/30 (*continued*)

Study	View	SID in/cm	kVp	mAs
Tarsus				
	DP	30/76	80	1.7
	Medial Oblique	30/76	80	1.7
Angle Down 5	Lateral	30/76	80	1.7
Angle Up 5	Lateral Oblique	30/76	80	1.7
If Indicated	Calcaneal Skyline	18/46	80	1.7
	Flexed DP	30/76	100	1.7
Skull (TMJ & Bullae)				
	DV	30/76	100	2.5
	Lateral	30/76	100	1.4
	Right Oblique	30/76	100	1.4
	Left Oblique	30/76	100	1.4
Sinus/Dental				
	DV	30/76	100	2.5
	Lateral	30/76	80	1.7
	Right Oblique	30/76	80	1.7
	Left Oblique	30/76	80	1.7
Cervical Spine				
	Cranial (Lateral)	30/76	100	1.7
	Middle (Lateral)	30/76	100	2
	Caudal (Lateral)	30/76	100	2.5
Pharynx	Lateral	30/76	80	1.7
Shoulder	Lateral	30/76	100	3
Elbow				
	Cran->Caud	30/76	100	1.4
	Medio-Lateral	30/76	100	1.7
Radius				
	Cran->Caud	30/76	80	1.7
	Lateral	30/76	80	1.7
	Obliques	30/76	80	1.7
Stifle				
	Caud->Cran	30/76	100	2
	Lateral	30/76	80	1.7
	Lateral Oblique	30/76	80	1.7
	Patellar Skyline	18/46	80	1.7
Tibia				
	Caud->Cran	30/76	100	1.4
	Lateral	30/76	80	1.7
	Obliques	30/76	80	1.7
Withers	Lateral	30/76	80	1.2

Appendix A: Digital Radiography Technique Charts *(continued)*

Gantry

Study	View	SID in/cm	kVp	mAs
Navicular				
	Navicular DP	40/102	65	5
	DP P3	40/102	80	2
	Foot DP	40/102	80	3.2
	Lateral	40/102	80	3.2
	Navicular Skyline	40/102	80	3.2
If Indicated	Nav. Obliques	40/102	80	2
If Indicated	P3 Obliques	40/102	80	2
P-3				
	Lateral	40/102	80	3.2
	DP P3	40/102	80	2
If Indicated	P3 Obliques	40/102	80	2
Pastern				
	DP	40/102	80	3.2
	Lateral	40/102	80	3.2
	Dlpmo	40/102	80	3.2
	Dmplo	40/102	80	3.2
Fetlock				
	DP	40/102	80	3.2
	Dlpmo	40/102	80	3.2
	Lateral	40/102	80	3.2
	Dmplo	40/102	80	3.2
	Flexed Lateral	40/102	80	3.2
If Indicated	Special Palmar	40/102	80	3.2
Metacarpus/Metatarsus				
	DP	40/102	80	3.2
	Lateral	40/102	80	3.2
	Medial Oblique	40/102	80	3.2
	Lateral Oblique	40/102	80	3.2
Carpus				
	DP	40/102	80	3.2
	Dlpmo	40/102	80	3.2
	Lateral	40/102	80	3.2
	Dmplo	40/102	80	3.2
	Flexed Lateral	40/102	80	3.2
If Indicated	Distal Rad. Sky	40/102	80	3.2
	Distal Row Sky	40/102	80	3.2
If Indicated	Prox. Row Sky	40/102	90	3.2

Appendix A: Digital Radiography Technique Charts

Gantry (continued)


Study	View	SID in/cm	kVp	mAs
Tarsus				
	DP	40/102	95	3.2
	Dlpmo	40/102	80	3.2
Angle Down 5	Lateral	40/102	80	3.2
Angle Up 5	Dmplo	40/102	80	3.2
If Indicated	Calcaneal Skyline	40/102	80	3.2
Skull (TMJ & Bullae)				
	DV	40/102	100	4
	Lateral	40/102	90	3.2
	Right Oblique	40/102	90	3.2
	Left Oblique	40/102	90	3.2
Sinus/Dental				
	DV	40/102	100	3.2
	Lateral	40/102	90	3.2
	Right Oblique	40/102	90	3.2
	Left Oblique	40/102	90	3.2
Cervical Spine				
	Cranial (Lateral)	40/102	90	3.2
	Middle (Lateral)	40/102	100	4
	Caudal (Lateral)	40/102	110	15
Pharynx	Lateral	40/102	90	3.2
Shoulder	Lateral	40/102	110	15
Elbow				
	Cran->Caud	40/102	110	4
	Medio-Lateral	40/102	90	3.2
Radius				
	Cran->Caud	40/102	100	3.2
	Lateral	40/102	90	3.2
Stifle				
	Caud->Cran	40/102	110	12
	Lateral	40/102	90	8
	Pldmo	40/102	95	3.2
Tibia				
	Caud->Cran	40/102	100	8
	Lateral	40/102	90	6
Withers	Lateral	40/102	100	4

Appendix B:

Mark X™ System Field Replacement Guide



1. Using a Phillips screwdriver remove four (4) screws.
2. Box will remain on because of 3 Guide Pins. Gently slide electronics box off Guide pins, By pulling back. There are no wires to disconnect.




Appendix B: Mark X™ System Field Replacement Guide (continued)


3. Flip electronics box over and pull out 16GB Compact Flash Card. The compact flash card will have the most recent Database file on it with all settings required to restore.



You will need to place the old Compact Flash into the new electronics box replacement.




To Slide out old 16GB CF via the cut out:
Gently pull left or right tab with thumb
while pulling card back towards "X".




4. You will need to place the old Compact Flash into the new electronics box replacement.

To Slide in new 16GB CF via the cut out:
Gently align card up with rails and push card
forwards towards "X".



Appendix B: Mark X™ System Field Replacement Guide (continued)

5. After replacing CF card into new Electronics Box slide on using guide pins.
The guide pins will line up the blind mate connector for you.



6. Using a Phillips screwdriver to tighten four (4) screws.



Appendix B: Mark X™ System Field Replacement Guide (continued)

7. Obtain Internet Connection, Power On, then call Technical Support at 1.800.819.5538 to restore settings. Restore procedure will take roughly 20 minutes.



Appendix C: Calibration Procedure

1. Launch the TruDR application as Administrator (username: admin password: admin).
2. Select the Actions button located in the menu bar at the top of the screen.
3. Select the **Manage** button located at the top of the dropdown that appears after selecting **Actions**.
4. Select **Equipment** from the menu on the left side of the screen.
5. The **Equipment Configuration** screen will appear in the window on the right. At the bottom of the window select the **Perform Calibration** button.



You may need to minimize the On-Screen Keyboard in order to see the button. To do this select the Keyboard button located in the top left corner of the screen.)

6. The **Calibration** screen will load on top of the application. To begin the calibration, select the **Perform Gain** button. Once you have done this, the application will begin to acquire a series of 4 **Dark** images 30 seconds apart and average them to create the **Offset** image. After each image is acquired, the Progress text box will read: Please Wait.
7. Next, you will be prompted to take a shot to check technique. The target value for the gain calibration is ~11000 ADU.



A good starting technique at 50" SID is 80 kVp and .17 seconds

8. If the setting is correct the **Calibration Progress** box will display the **Average Exposure Count** and allow you to continue the calibration. If the settings are not correct, the **Calibration Progress** box will tell you if the settings are too high or too low, reset, and then prompt you to adjust your settings and take another image. (The window for acceptable values is +/- 20%.)
9. If your settings were correct, you will be prompted to acquire a series of 4 images 30 seconds apart.
10. Once all 4 gain calibration images have been created, you will be prompted to acquire a **Flat Field** image. Acquire this image at the same technique as your calibration images.
11. Once the **Flat Field** image has been acquired, the application will begin calculating the statistics of the calibration. **DO NOT CLICK ANY BUTTONS OR EXIT THE SCREEN UNTIL THE STATISTICS HAVE BEEN GENERATED.**

Appendix C: Calibration Procedure *(continued)*

12. Approximate Calibration Statistics:
 - Offset Median: 2129
 - Offset Std Dev: 258.1529
 - Gain Median: 9586
 - Gain Std Dev: 144.6174
13. Using the dropdown window next to the **Get Image** button at the bottom of the Calibration Screen you can view the **Offset**, **Gain**, and **Flat-Field** images by selecting them in the dropdown and selecting the **Get Image** button. Use the **Brightness** and **Contrast** sliders at the bottom of the image window to adjust the image.
14. Select the **Flat-Field** image from the dropdown window and select **Get Image**. Inspect the image for any defective lines, pixels, or any other artifact. If the image is a mottled/uniform gray color then the calibration was successful.

Appendix D: Technical Support

Shipping Address

Sound-Eklin™

6359 Paseo Del Lago

Carlsbad . California 92011

USA

Technical Support

Office hours: Weekdays 8:00 A.M. — 5:00 P.M. Pacific Time

Emergency 24-hour support is available.

Telephone: 800-268-5354 option 3

Website

www.soundeklin.com



a VCA ANTECH company

Sound-Eklin™ Technical Support

6359 Paseo Del Lago
Carlsbad . California 92011 USA

toll free: 800.268.5354 option 3
phone: 760.918.9626
fax: 76.918.9620
international: +1.760.918.9626

www.soundeklin.com